



heinlein
PLASTIK-TECHNIK



flipdropper[®]

Flip top cap with integrated drip function

||||| A MEMBER OF ALPLA INDEN PHARMA |||||



OUR PRODUCT INNOVATION 2020

SOPHISTICATED DOWN TO THE LAST DETAIL



Visible Tamper-evident closure



No unscrewing necessary



Tamper-evident seal remains in the closure



Perfect dosage



flipdropper®

NEW GENERATION OF PHARMA PACKAGING

As an innovative manufacturer of reliable primary packaging materials, we have developed a new product:

**The flipdropper® –
single-piece flip top cap with tamper-evidence
and integrated drip function**

The special feature:

The flip top cap can be opened in one easy step after which the contents can be dispensed. Cumbersome unscrewing is now a thing of the past. This is how the Heinlein flipdropper® combines user-friendly handling with ideal product protection.

STEP 1:

PUSH UP THE TOP



Simply press upwards to open

Upon initial opening, the seal falls into a pocket provided for this purpose. No unscrewing necessary.

STEP 2:

TURN DOWNWARDS & DROP



Intuitive operation and perfect dosage

The integrated universal drip function starts dripping as soon as the container is turned down.

STEP 3:

PULL DOWN THE TOP



Safe sealing

Can be securely resealed thanks to the intuitive push/pull mechanism.



flipdropper®

SHORT FACTS

Benefits

- + no tiresome screwing on and off
- + quick product access, simple push/pull mechanism
- + perfect dosage
- + for liquids with different viscosities
- + tamper-evident closure - ideal product protection
- + can be processed on all commonly encountered sealing machines
- + Pharma quality: Material PE-HD

Application areas: Pharmacy and medicine, food supplements, lifestyle products

Item number: KVT-18.202

Thread sizes: DIN 168 / GL 18

Award winner:



Please contact us for samples.



WE PROTECT YOUR VALUES

PHARMA PACKAGING MADE IN GERMANY